



Recombinant Mouse WD repeat-containing protein 24 (Wdr24)

Product Code	CSB-MP804957MO
Abbreviation	Wdr24
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CFJ9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MEKMSRVSTA LSGSALTGRT MHCHLDAPAN AISVCRDAAQ VVAGRSIFK IYAIEEEQFV EKLNLRVGRK PSLNLSCADV VWHQMDENLL ATAATNGVVV TWNLGRPSRN KQDQLFTEHK RTVNKVCFHP TEAHVLLSGS QDGFMKCFDL RRKDSVSTFS GQSESVRDVQ FSIRDYFTFA STFENGNVQL WDIRRPDRCE RMFTAHNQPV FCCDWHPEDR GWLATGGRDK MVKVWDMTTH RAKEIHCVQT IASVARVKWR PECRHHLATC SMMVDHNIYV WDVRRPFVPA AMFEEHRDVT TGIAWRHPHD PSFLLSGSKD STLCQHLFRD ASQPVERANP EGLCYGLFGD LAFAVKESLV AESGRKPYA GDRRHPIFFK RKLDP AEPFS GLASSALSVF ETESSGGSMS WFDVTAERYV LAGRPLAELC DHNAKVAREL GRNQVAQTWT MLRIIYCS PG LVSSANLNHS VGKGSSCGLP LMNSFNLKDM GPGLGSETRL DRSKGDTRSD AALLDSSATL VTNEDNEETE GSDVPADYLL GDVEGEDDEL YPLDTEHVHS EEPEYVLPQE AFPLRHEIVD TSPGPEHLQD KADSPHVSGN EADTASLAPV DSSSSLLSVS HALYDSRLPP DFFSVLVDRM LRFYAEQGDV QMAVSVLIVL GERVRKDIDE QTQEHWYTSY IDLLQRFC LW NVSNEVVKLS TSRAVSCLNQ ASTTLHVNC S HCKRPMSSRG WVCDRCHRCA SMCVCHHV KGLFVWCQGC SHGGHLQHIM KWLEGSSHCP AGCGHLCEYS
Source	Mammalian cell
Target Names	Wdr24
Protein Names	Recommended name: WD repeat-containing protein 24
Expression Region	1-790
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.