



# Recombinant Mouse WD repeat-containing protein 24 (Wdr24)

<b>Product Code</b>	CSB-YP804957MO
<b>Abbreviation</b>	Wdr24
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8CFJ9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEKMSRVSTA LSGSALTGRT MHCHLDAPAN AISVCRDAAQ VVAGRSIFK IYAIEEEQFV EKLNLRVGRK PSLNLSCADV VWHQMDENLL ATAATNGVVV TWNLGRPSRN KQDQLFTEHK RTVNKVCFHP TEAHVLLSGS QDGFMKCFDL RRKDSVSTFS GQSESVRDVQ FSIRDYFTFA STFENGNVQL WDIRRPDRCE RMFTAHNQPV FCCDWHPEDR GWLATGGRDK MVKVWDMTTH RAKEIHCVQT IASVARVKWR PECRHHLATC SMMVDHNIYV WDVRRPFVPA AMFEEHRDVT TGIAWRHPHD PSFLLSGSKD STLCQHLFRD ASQPVERANP EGLCYGLFGD LAFVAVKESLV AESGRKPYA GDRRHPIFFK RKLDPAPFP GLASSALSVF ETESSGGSMS WFDVTAERYV LAGRPLAELC DHNAKVAREL GRNQVAQTWT MLRIIYCS PG LVSSANLNHS VGKGSSCGLP LMNSFNLKDM GPGLGSETRL DRSKGDTRSD AALLDSSATL VTNEDNEETE GSDVPADYLL GDVEGEDDEL YPLDTEHVHS EEPEYVLPQE AFPLRHEIVD TSPGPEHLQD KADSPHVSGN EADTASLAPV DSSSSLLSVS HALYDSRLPP DFFSVLVDRM LRFYAEQGDV QMAVSVLIVL GERVRKDIDE QTQEHWYTSY IDLLQRFCW NVSNEVVKLS TSRAVSLNQ ASTTLHVNC S HCKRPMSSRG WVCDRCHRCA SMCVAVCHHV KGLFVWCQGC SHGGHLQHIM KWLEGSSHCP AGCGHLCEYS
<b>Source</b>	Yeast
<b>Target Names</b>	Wdr24
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 24
<b>Expression Region</b>	1-790
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.