



# Recombinant Buchnera aphidicola subsp. Baizongia pistaciae Transcription termination factor Rho (rho)

<b>Product Code</b>	CSB-EP805054BMZ
<b>Abbreviation</b>	rho
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q89A22
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Buchnera aphidicola subsp. Baizongia pistaciae (strain Bp)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNLTALKNIP VSELIFLGDN AGLLENLARMR KQDIIFSILK QHAKSGEDIF GDGVLEILQD GFGFLRSSDS SYLAGPDDIY VSPSQIRRFN LRTGDTISGK IRPPKEGERY FALLKVNKVN YDKPENARSK ILFENLTPH ANSRLRMERG NGSTEDLTAR VLDLAPIGR GQRGLIVAPP KAGKTMLLQN IAQSIAYNHP DCVLMVLLID ERPEEVTEMR RLVKGEVVAS TFDEPSSRHV QVSEMVIEKA KRLVEHKKDV IILLDSITRL ARAYNTVPA SGKVLTGGVD ANALHRPKRF FGAARNVEEG GSLTIIATAL IDTGSKMDEV IYEEFKGTGN MELPLSRKIA EKRVFPAIDY NRSCTRKEEL LTKSDELQKM WILRKIIHPM SEIDAMEFLI NKLAMTKTNN EFFDMMKRS
<b>Source</b>	E.coli
<b>Target Names</b>	rho
<b>Protein Names</b>	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
<b>Expression Region</b>	1-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.