



Recombinant Danio rerio RuvB-like 1 (ruvbl1)

Product Code	CSB-BP805217DIL
Abbreviation	ruvbl1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8AWW7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKIEEVKSTT KTQRIASHSH VKGLGLDEAG NAKQSASGLV GQESAREACG IITELIRSKK MAGRAILLAG PPGTGKTALA LAMAQELGNK VPFCPMVGSE VYSSEIKKTE VLMENFRAI GLRIKETKEV YEGETELTP CETENPMGGY GKTISHVIIG LKTAKGTKQL KLDPSIYESL QKERVEVDV IYIEANSGAV KRQGRCDTFA TEFDLAEEY VPLPKGDVHE KKEIIQDVTL HDLDVANARP QGGQDILSMM GQLMKPKKTE ITDKLRGEIN KVVNKYIDQG VAELVPGVLF IDEVHMLDIE CFTYLHRALE SSIPIVVFA SNRGNCLIRG TEDISSPHGI PLDLLDRVMI IRTMLYTPQE MKQIIKIRAQ TEGLNISEEA LSHLGEIGTK TTLRYAVQLL TPASLLARVQ GREVVEKEHV EEINELFYDA KSSAKILQDQ HTKFMK
Source	Baculovirus
Target Names	ruvbl1
Protein Names	Recommended name: RuvB-like 1 EC= 3.6.4.12 Alternative name(s): Pontin zPontin
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.