



# Recombinant Mouse Pre-mRNA-splicing factor SLU7 (Slu7)

<b>Product Code</b>	CSB-EP805272MO
<b>Abbreviation</b>	Slu7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BHJ9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SAAAVDPVS ATPMTGSKEM SLEEPKKMTR EDWRKKKELE EQRKLGNAPA EVDEEGKDIN PHIPQYISSV PWYIDPSKRP TLKHQRPQPE KQKQFSSSGE WYKRGVKENS ITTKYRKGAC ENCGAMTHKR KDCFERRRV GAKFTGTNIA PDEHVQPQLM FDYDGKRDRW NGYNPEEHMK IVEEYAKVDL AKRTLKAQKL QEELASGKLV EQANSPKHQW GEEEPNSQME KDHNSEDEDE DKYADDIDMP GQNFDSKRRI TVRNLRNIRD IAKYLRNLDP NSAYYDPKTR AMRENPYANA GKNPDEVSYA GDNFVRYTGD TISMAQTQLF AWEAYDKGSE VHLQADPTKL ELLYKSFVKV KEDFKEQQKE SILEKYGGQE HLDAPPAELL LAQTEDYVEY SRHGTVIKGG ERAVACSKYE EDVKINNHTH IWGSYWKEGR WGYKCCHSFF KYSYCTGEAG KESVNSEECI ITGATAEESV KKPQALLELH QEKLKEEKKK KKKKKKHRKS SSSDDEERK QEKLKALNA EEARLLHVKE IMQIDERKRP YNSIYETREP TEEEMEAYRM KRQRPDDPMA SFLGQ
<b>Source</b>	E.coli
<b>Target Names</b>	Slu7
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor SLU7
<b>Expression Region</b>	2-585
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Pre-mRNA splicing occurs in two sequential transesterification steps. This protein is a splicing factor that has been found to be essential during the second catalytic step in the pre-mRNA splicing process. It associates with the spliceosome and contains a zinc knuckle motif that is found in other splicing factors and is involved in protein-nucleic acid and protein-protein interactions.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.