



Recombinant Mouse Casein kinase I isoform alpha (Csnk1a1)

Product Code	CSB-BP805291MO
Abbreviation	Csnk1a1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BK63
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ASSSGSKAE FIVGGKYKLV RKIGSGSFGD IYLAINITNG EEVAVKLESQ KARHPQLLYE SKLYKILQGG VGIPHIRWYG QEKDYNVLVM DLLGPSLEDL FNFCRRFTM KTVLMLADQM ISRIEYVHTK NFIHRDIKPD NFLMGIGRHC NKLFLIDFGL AKKYRDNRTR QHIPYREDKN LTGTARYASI NAHLGIEQSR RDDMESLGYV LMYFNRTSLP WQGLKAATKK QKYEKISEKK MSTPVEVLCK GFPAEFAMYL NYCRLRFEE APDYMRLRQL FRILFRTLNH QYDYTFDWTM LKQKAAQQA SSSGQQQAQ TPTGKQTDKT KSNMKGK
Source	Baculovirus
Target Names	Csnk1a1
Protein Names	Recommended name: Casein kinase I isoform alpha Short name= CKI-alpha EC= 2.7.11.1 Alternative name(s): CK1
Expression Region	2-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.