



# Recombinant Mouse Brain-specific angiogenesis inhibitor 1-associated protein 2 (Baiap2)

<b>Product Code</b>	CSB-MP805295MO
<b>Abbreviation</b>	Baiap2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BKX1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLSLRSEEMH RLTENYKTI MEQFNPSLRN FIAMGKNYEK ALAGVTFAAK GYFDALVKMG ELASESQGSK ELGDVLFQMA EVHRQIQNQL EETLKSFHNE LLTQLEQKVE LDSRYLSAAL KKYQTEQRSK GDALDKCQAE LKCLRKKSQG SKNPQKYSK ELQYIDAISN KQGELENYVS DGYKTALTEE RRRFCFLVEK QCAVAKNSAA YHSGKELLA QKLPLWQQAC ADPNKIPDRA VQLMQQMANS NGSILPSALS ASKSNLVID PIPGAKPLPV PPELAPFVGR MSAQENVPVM NGVAGPDSYD YNPWADRKAA QPKSLSPQSQ QSKLSDSYSN TLPVRKSVTP KNSYATTENK TLPRSSSMAA GLERNRMRV KAIFSHAAGD NSTLLSFKEG DLITLLVPEA RDGWHYGESE KTKMRGWFPF SYTRVLDSYG SDRLHMSLQQ GKSSSTGNLL DKDDLALPPP DYGTSSRAFP TQTAGTFKQR PYSVAVPAFS QGLDDYGARS VSRNPFANVH LKPTVTNDRS APLLS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Baiap2
<b>Protein Names</b>	Recommended name: Brain-specific angiogenesis inhibitor 1-associated protein 2 Short name= BAI-associated protein 2 Short name= BAI1-associated protein 2 Alternative name(s): Insulin receptor substrate protein of 53 kDa Short nam
<b>Expression Region</b>	1-535
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.