



Recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 4 (Adamts4)

Product Code	CSB-EP805312MO-B
Abbreviation	Adamts4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BNJ2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	FA SLSRFVETLV VADDKMAAFH GTGLKRYLLT VMAAAAKAFK HPSIRNPVNL VVTRLVILGS GQEGPQVGPS AAQTLRSFCT WQRGLNTPND SDPDHFDTAI LFTRQDLGCV STCDTLGMAD VGTVCDPARS CAIVEDDGLQ SAFTAAHELHGVFVFNMLHDNS KPCTNLNGQG GSSRHVMAPV MAHVDPEEPW SPCSARFITD FLDNGYGHCL LDKPEAPLHL PATFPGKDYD ADRQCQLTFG PDSSHCPQLP PPCAALWCSG HLNGHAMCQT KHSPWADGTP CGSSQACMGG RCLHVDQLKD FNVPQAGGWG PWGPWGDCSR TCGGGVQFSS RDCTRPVPRN GGKYCEGRRT RFRSCNTENC PHGSALTFRE EQCAAYNHRT DLFKSFPGPM DWVPRYTGVA PRDQCKLTCQ ARALGYYYVL EPRVADGTPC SPDTSSVCVQ GRCIHAGCDR IIGSKKKFDK CMVCGGDGSR CSKQSGSFKK FRYGYSDVVT IPAGATHILV RQQGGGSLKS IYLALKLSDG SYALNGEYTL MPSPDVTVLP GAVSLRYSGA TAASETLSGH GPLAQPLTLQ VLVAGNPQNA RLRYSFFVPR PVPSTPRPPP QDWLQRRAEI LKILRKRPA GRK
Source	E.coli
Target Names	Adamts4
Protein Names	Recommended name: A disintegrin and metalloproteinase with thrombospondin motifs 4 Short name= ADAM-TS 4 Short name= ADAM-TS4 Short name= ADAMTS-4 EC= 3.4.24.82 Alternative name(s): Aggrecanase-1
Expression Region	209-833
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the ADAMTS (a disintegrin and



metalloproteinase with thrombospondin motifs) protein family. Members of the family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The enzyme encoded by this gene lacks a C-terminal TS motif. It is responsible for the degradation of aggrecan, a major proteoglycan of cartilage, and brevican, a brain-specific extracellular matrix protein. The cleavage of aggrecan and brevican suggests key roles of this enzyme in arthritic disease and in the central nervous system, potentially, in the progression of glioma.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.