



Recombinant Mouse Serine/threonine-protein kinase PAK 7 (Pak7)

Product Code	CSB-MP805374MO
Abbreviation	Pak7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8C015
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MFGKKKKKIE ISGPSNFEHR VHTGFDPQEQ KFTGLPQQWH SLLADTANRP KPMVDPSCIT PIQLAPMKTI VRGNKCKET SINGLLEDFD NISVTRSNSL RKESPPTPDQ GAASRIQGH S EENGFITFSQ YSSESDTTAD YTTEKYRDRS LYGDDL DLYY KSSHA AKQNG HAMKMKHGDA YPEMKS LKT DLAGFPVDYH THLDSL RKSS EYGDLRWDYQ RASSSSPLDY SFQLTPSR TA GTSRCSKESL AYSESDWGPS LDDYDRR PKS SYLHQTSPQP AMRQRSKSGS GLQEPMPFG ASAFKTHPQG HSYNSYTYPR LSEPTMCIPK VD YDRAQMV F SPPLSGSDTY PRGPTKLPQS QSKAGYSSGS HQYPSGYHKA SLYHHPSLQT SSQYISTASY LSSLSISSST YPPPSWGSSS DQQPSRV SHE QFRAALQLVV SPGDPREYLD NFIKIGEGST GIVCIATEKH TGKQVAVK KM DLRKQQRREL LFNEVVIMRD YHHDNVDMY NSYLVGDELW VVMEFLE GGA LTDIVTHTRM NEEQIATVCL SVLKALSYLH NQGV IHRDIK SDSILLTSDG RIKLSDFGFC AQVSKEVPKR KSLVGTPYWM APEVISRLPY GTEVDIWSLG IMVIEMIDGE PPYFNEPPLQ AMRRIRD SLP PRVKDLHKVS SMLRGFLDLM LVREPSQRAT AQELLGHPFL KLAGPPSCIV PLMRQYRHH
Source	Mammalian cell
Target Names	Pak5
Protein Names	Recommended name: Serine/threonine-protein kinase PAK 7 EC= 2.7.11.1 Alternative name(s): p21-activated kinase 5 Short name= PAK-5 p21-activated kinase 7 Short name= PAK-7
Expression Region	1-719
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the PAK family of Ser/Thr protein kinases. PAK family members are known to be effectors of Rac/Cdc42 GTPases, which have



been implicated in the regulation of cytoskeletal dynamics, proliferation, and cell survival signaling. This kinase contains a CDC42/Rac1 interactive binding (CRIB) motif, and has been shown to bind CDC42 in the presence of GTP. This kinase is predominantly expressed in brain. It is capable of promoting neurite outgrowth, and thus may play a role in neurite development. This kinase is associated with microtubule networks and induces microtubule stabilization. The subcellular localization of this kinase is tightly regulated during cell cycle progression. Alternatively spliced transcript variants encoding the same protein have been described.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.