



Recombinant Mouse Smad nuclear-interacting protein 1 (Snip1)

Product Code	CSB-BP805905MO
Abbreviation	Snip1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BIZ6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MKAGKSERER SGRRRHRSGD ALTTVVVKQE RLSPEPVAHR RPDAPAASLS PPAAEPGHSG HRGSRARSPA KKKSKSSGRR SKSPRTKRSQ SPHYPMVKVK QEREDHPRRG REDRQHREPS EQEHRRARNS ERDRHRGHSR QGRSSDERPV SGQDRDRDSQ NLQAQEEERD FHNARRREHR QQNESAGSEA QEVIPRPAGN RSKEVPVKEK PSFELSGALL EDTNTFRGVV IKYSEPPEAR IPKKRWRLYP FKNDEVLPVM YIHRQSAYLL GRHRIADIP IDHPSCSKQH AVFQYRLVEY TRADGTVGRR VKPYIIDLGS GNGTFLNNKR IEPQRYELK EKDVLFKFGFS SREYVLLHES SDTSELDRKE DEDDEEEEMV SDS
Source	Baculovirus
Target Names	Snip1
Protein Names	Recommended name: Smad nuclear-interacting protein 1
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.