



Recombinant Mouse Interleukin-33 (Il33), partial

Product Code	CSB-EP805971MO-B
Abbreviation	Il33
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8BVZ5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SI QGTSLLTQSP ASLSTYNDQS VSFVLENGCY VINVDDSGKD QEQQVLLRY YESPCPASQS GDGVDGKKLM VNMSPIKDTD IWLHANDKDY SVELQRGDVS PPEQAFFVLH KKSSDFVSFE CKNLPGTYIG VKDNQLALVE EKDESCNNIM FKLSKI
Source	E.coli
Target Names	Il33
Protein Names	Recommended name: Interleukin-33 Short name= IL-33
Expression Region	109-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>