



# Recombinant Mouse FERM domain-containing protein 6 (Frmd6)

<b>Product Code</b>	CSB-YP806006MO
<b>Abbreviation</b>	Frmd6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8C0V9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNKLTFHNNK AMQDRRRVCI FLPNDKSVSI IINVKILCHQ LLVQVCDLLR LKDSHLFGLS VIQNEHVYM ELSQKLYKYC PKEWKKEASK VRQYEVTWGI DQFGPPMIH FRVQYYVENG KLISDRIARY YYYWHLRKQV LHSQCVLREE AYFLLAAFAL QADLGNFKRK LHHGDYFEPE AYFPAWVVS KRGKDYILKHI PNMHKDQFAL TASEAYLKYI KEAVRLDDVA IHYYRLYKDK REAEGSLTLG LTMRGIQIFQ NLEEEKQLLY DFPWTVNGKL VFGVGGKFEIL PDGLPSARKL VYYTGCPTRS RHLLQLLSNS HRLYMNLPV LRHLRKQEEEN EEKKQYRESY ISDNLDLMD PLEKRSRSG SSAGSVKHKR LSRHSTASHS SSHTSGIEAD TKPRDPGPEP SCSSGSAMHRK LKTCSSMTSH GSSHTSGVES GGKDRLEEDS QDEEIEMLVD DPRDLEPMPE ESLEVSPMC IYITEDMLLS RKLNGHSGLI VKEIGSSTSS SSETVVRLRG QSTDSPQTI CRKPKTSTDR HLSLDDIRL YQKDFLRIAG LCQDTAQS YT FCGGHELDES GLYCNSCLAQ QCVNIQDAFP VKRASKYFSL DLTHDEVPEF VV
<b>Source</b>	Yeast
<b>Target Names</b>	Frmd6
<b>Protein Names</b>	Recommended name: FERM domain-containing protein 6
<b>Expression Region</b>	1-622
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.