



Recombinant Mouse Coiled-coil domain-containing protein 63 (Ccdc63)

Product Code	CSB-EP806050MO-B
Abbreviation	Ccdc63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CDV6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MKKVRRKGSS GLSESEKAR EQLAQAE LRK LRQQFRKMVD SRKSFNFRHQ RMIAGQYKEI ETLKAEQAET TMLLSLVKSP KNLDINQKNF MELRFLQTK GDYEALISSM KVLLGELDDK IVQMERKITN QRQLFLRTQE ANNPKRLQRQ IHVLESRLNL VTVHFDTMLT SNSQLRKEIE DLLFEKAAYD HVYQQLQRRL QTQKKT MNVA IEQSAQAYEQ RVEAMARMAA MKDRQQKDIS QYNLEIRELE RLYDHETK LK SFLAKLNDR SEFEDQAKKQ EDVKS KLGK KGKGFAS Y EVAHLRLLKL AENGLNQLT EDFLAKEEKN FARFTYVTEL NNDMETMHKK TQRIQDDIIN LRSQQQTSHE GTRSILKQME EKLRKTTQEA DIYETKYKEM SKTLEYLKNS VEKMFKKINC DATEILGKLG ESGKVTDINL QQYFAIIEKK TNDLLLLESF RRLQEAEGPD VDVPPFVNP FWGGSALLKP PEPIRVVLPV FGADSFSDKL EEVDSPLDHS TLQQMVLENF LQRERTKELQ DTMSEK GDEI RLKKKVIG
Source	E.coli
Target Names	Ccdc63
Protein Names	Recommended name: Coiled-coil domain-containing protein 63
Expression Region	1-558
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.