



Recombinant Rat Thiamine-triphosphatase (Thtpa)

Product Code	CSB-YP806072RA
Abbreviation	Thtpa
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CGV7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	AQGLIEVER KFTPGPDTEE RLQKLGATLE HRVTFRDITYY DTSELSLMLS DHWLRQREGS GWEFKCPGVT GVSOPHNEYV EVTSESAIVT QLFELLGSGE QETAGVA AVL GRLKLQEVAS FITTRSSWKL ALSGAHEEES LLTVDL DSTD FGYAVGEVEA VVHEKAEVPA ALEKIISVSS MLGVPAQEKA PAKLLVYLQR FRPQDYQRLL EADSSGEATG DSV P
Source	Yeast
Target Names	Thtpa
Protein Names	Recommended name: Thiamine-triphosphatase Short name= ThTPase EC= 3.6.1.28
Expression Region	2-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.