



# Recombinant Chicken Endophilin-A2 (SH3GL1)

<b>Product Code</b>	CSB-YP806325CH
<b>Abbreviation</b>	SH3GL1
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q8AXV0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSVAGLKKQF YKASQLVSEK VGGAEGTKLD DDFKEMEKKV DLTSKAVTEV LTRTIEYLQP NPASRAKLTM LNTMSKIRGQ VKNPGYPQSE GLLGESMIRY GKELGEDSNF GDALLDAGES MKRLAEVKDS LDIEVKQNF I DPLQNLCDKD LKEIQHHLKK LEGRRLDFDY KKKRQGKIPD EELRQAMEKF EESKEVAETS MHNLLET DIE QVSQLSALVD AQLDYHRQAV QILDELAEKL KRRMREASSR PRREYKPKPR ETYDFGESDQ SNGGF SCTPT PKVSASSSFR SDKPFRTSVR SIPHLDQPCC KALYDFEPEN DGELGFKEGD IITLTNQIDE NWEYEGMINGQ SGFFPLNYVE VLVPLPQ
<b>Source</b>	Yeast
<b>Target Names</b>	SH3GL1
<b>Protein Names</b>	Recommended name: Endophilin-A2 Alternative name(s): Endophilin-2 SH3 domain-containing GRB2-like protein 2 SH3p8
<b>Expression Region</b>	1-367
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>