



Recombinant Vesicular stomatitis Indiana virus

Matrix protein (M)

Product Code	CSB-EP806800VBF-B
Abbreviation	M
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8B0I2
Product Type	Recombinant Protein
Immunogen Species	Vesicular stomatitis Indiana virus (strain 98COE North America) (VSIV)
Purity	≥85% (SDS-PAGE)
Sequence	MSSLKKILGL KGKGKKSKKL GIAPPPYEED TSMEYAPSAP IDKSYFGVDE MDTHDPNQLR YEKFFFTVKL TVRSNRPFRT YSDVAAAVSH WDHMYIGMAG KRPFYKILAF LGSSNLKATP AVLADQGQPE YHAHCEGRAY LPHRMGKTPP MLNVPEHFRR PFNIGLYKGT IELTMTIYDD ESLEAAPMIW DHFNSSKFSD FREKALMFGL IVEKKASGAW ILDSVSHFK
Source	E.coli
Target Names	M
Protein Names	Recommended name: Matrix protein
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.