



Recombinant Mouse Spermatogenesis-associated protein 2-like protein (Spata2l)

Product Code	CSB-MP806882MO
Abbreviation	Spata2l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BNN1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MGSSSLSEDY RQCLERELRR GRAGVCGDPS LRAVLWQILV EDFDLHGALQ DDALALFTDG LWGRADLAPA LQDLARAFEL LELAAVHLYL LPWRKEFTTI KTFSGGYVHV LKGVLSEELL TRSFQKMGYV RRDNHRLMVT TPPPACQLVQ VALGCFALRL ECEILSEVLT QLGTSVLP AE ELLRARRASG DVASCVAWLQ QRLAQDEEPP PLPPRGTPAT YGAPVDLYQD LQEDESSEAS LYGEPSPGLD SPPVELAYRP PLWEQSAKLW GSGGQPWEPP ADDMHRASSP PYGALEEELE PEPSAFSFLS LRRELSRSGD LAPPESPSSP GQASPRHRQA EAAASSAYGP AVEPLSYQAH SCLSPGNLPT LCCDTCRQLH ATHCTALSAC RPTHSLRILL GDNQRRLWLQ RAQVDNLLYD SPGAHP
Source	Mammalian cell
Target Names	Spata2l
Protein Names	Recommended name: Spermatogenesis-associated protein 2-like protein Short name= SPATA2-like protein
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.