



Recombinant Mouse Spermatogenesis-associated protein 2-like protein (Spata2l)

| | |
|--------------------------|--|
| Product Code | CSB-MP806882MO |
| Abbreviation | Spata2l |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8BNN1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MGSSSLSEDY RQCLERELRR GRAGVCGDPS LRAVLWQILV EDFDLHGALQ DDALALFTDG LWGRADLAPA LQDLARAFEL LELAAVHLYL LPWRKEFTTI KTFSGGYVHV LKGVLSEELL TRSFQKMGYV RRDNHRMLMT TPPACQLVQ VALGCFALRL ECEILSEVLT QLGTSVLP AE ELLRARRASG DVASCVAWLQ QRLAQDEEPP PLPPRGTPAT YGAPVDLYQD LQEDESSEAS LYGEPSPGLD SPPVELAYRP PLWEQSAKLW GSGGQPWEPP ADDMHRASSP PYGALEEELE PEPSAFSFLS LRRELSRSGD LAPPESPSSP GQASPRHRQA EAAASSAYGP AVEPLSYQAH SCLSPGNLPT LCCDTCRQLH ATHCTALSAC RPTHSLRILL GDNQRRLWLQ RAQVDNLLYD SPGAHP |
| Source | Mammalian cell |
| Target Names | Spata2l |
| Protein Names | Recommended name: Spermatogenesis-associated protein 2-like protein Short name= SPATA2-like protein |
| Expression Region | 1-426 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.