



Recombinant Mouse TBC1 domain family member 17 (Tbc1d17)

Product Code	CSB-EP806934MO
Abbreviation	Tbc1d17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BYH7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MEGSSYRVVF EKGGVYLHTS ARKHQDPDSL IAGVIRVVEK DSDVFLHWAP VEEAGDPTQI LFKKDPSRGE PSTSEEEPTF DPGYEPDWAV ISTVRPQPHL AEPRKGAEPS SSRSSWAFSV SLGELKSIRR SKPGLSWAYL VLVTQAGGSL PALHFHRGGT RALLRVLSRY LLLASSPQDS RLYLVFPQDP SALSDSFHHL QLFDQDSSNV VSRFLQDPYS TTFSSFSRVT NFFRGALQPH PEGASSPNLP PLPDDEPEPG FEVISCVELG QRPTVERAPP VTEEEWNRVYV GPEGRLQNVP ELKNRIFSGG LSPGLRREAW KFLGGLYSWE SSAEEHKAHV RKKTDHEYFRM KLQWKSVAE QERRNSLLHG YRSLIERDVS RTDRTNKFYE GPENPGLSLL HDILLTYCMY HFDLGYVQGM SDLLSPILFV VQNEVDAFWC FCGFMELVHG NFEESQETMK RQLGQLLLLL RVLQDPLCDF LDSQDSGSLC FCFRWLLIWF KREFPFPDVL RLWEVLWTGL PGPNLHLLVA CAILDMERDT LMLSGFGSNE ILKHINELTM KLSVEDVLTR AEALYRQLTA CPELPHNVQE ILGLAQPEEP SSPSPVSPM PLSPTRAPLP PPLPEEVIPQ PDSSLEILPE DEDGA
Source	E.coli
Target Names	Tbc1d17
Protein Names	Recommended name: TBC1 domain family member 17
Expression Region	1-645
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.