



# Recombinant Mouse TBC1 domain family member 17 (Tbc1d17)

<b>Product Code</b>	CSB-MP806934MO
<b>Abbreviation</b>	Tbc1d17
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BYH7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEGSSYRVVF EKGGVYLHTS ARKHQDPDSL IAGVIRVVEK DSDVFLHWAP VEEAGDPTQI LFKKDPSRGE PSTSEEEPTF DPGYEPDWAV ISTVRPQPHL AEPRKGAEPS SSRSSWAFSV SLGELKSIRR SKPGLSWAYL VLVTQAGGSL PALHFHRGGT RALLRVLSRY LLLASSPQDS RLYLVFPQDP SALSDSFHHL QLFDQDSSNV VSRFLQDPYS TTFSSFSRVT NFFRGALQPH PEGASSPNLP PLPDDEPEPG FEVISCVELG QRPTVERAPP VTEEEWNRVYV GPEGRLQNVP ELKNRIFSGG LSPGLRREAW KFLGGLYSWE SSAEEHKAHV RKKTDHEYFRM KLQWKSVAE QERRNSLLHG YRSLIERDVS RTDRTNKFYE GPENPGLSLL HDILLTYCMY HFDLGYVQGM SDLLSPILFV VQNEVDAFWC FCGFMELVHG NFEESQETMK RQLGQLLLLL RVLQDPLCDF LDSQDSGSLC FCFRWLLIWF KREFPFPDVL RLWEVLWTGL PGPNLHLLVA CAILDMERDT LMLSGFGSNE ILKHINELTM KLSVEDVLTR AEALYRQLTA CPELPHNVQE ILGLAQPEEP SSPSPVSPM PLSPTRAPLP PPLPEEVIPQ PDSSLEILPE DEDGA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Tbc1d17
<b>Protein Names</b>	Recommended name: TBC1 domain family member 17
<b>Expression Region</b>	1-645
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.