



Recombinant Mouse Beclin 1-associated autophagy-related key regulator (Atg14)

Product Code	CSB-BP807000MO
Abbreviation	Atg14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CDJ3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MASPSGKGSW TPEAPGFGPR ALARDLVDSV DDAEGLYVAV ERCPLCNTTR RRLTCAKCVQ SGDFVYFDGR DRERFIDKKE RLSQLKNKQE EFQKEVLKAM EGKRLTDQLR WKIMSCKMRI EQLKQTICKG NEEMKKNSEG LLKNKEKNQK LYSRAQRHQE KKEKIQRHNR KLGDLVEKKT IDLKSHYERL ARLRRSHILE LTSIIFPIDE VKTSGRDPAD VSSETDSAMT SSMVSKLAEA RRTTYLSGRW VCDDHNGDTS ISITGPWISL PNNGDYSAYY NWVEEKTTQ GPDMEHNNPA YTISAALGYA TQLVNIVSHI LDINLPKCLC NSEFCGENLS KQKLTRAVRK LNANILYLC S QHVNLDQLQ PLHTLRNLMH LVSPRSEHLG RSGPFVEVRAD LEESMEFVDP GVAGESDASG DERVSDEETD LGTDWENLPS PRFCDIPSQP VEVSQSQSTQ VSPPIASSSA GGMISAAAS VTSWFKAYTG HR
Source	Baculovirus
Target Names	Atg14
Protein Names	Recommended name: Beclin 1-associated autophagy-related key regulator Short name= Barkor Alternative name(s): Autophagy-related protein 14-like protein Short name= Atg14L
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.