



Recombinant Mouse Actin-binding LIM protein 2 (Ablim2)

Product Code	CSB-EP807374MO-B
Abbreviation	Ablim2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BL65
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSAVSQPQAA HAPLEKPAST AILCNTCGNV CKGEVLRVQN KYFHIRCFVC KACGCDLAEG GFFVRQGEHI CTRDYQRLYG TRCFSCDRFI EGEVVSALGK TYHPDCFVCA VCRLPFPPGD RVTFNGKECM CQKCSPTLL GNSAHVAQGL RSCGGCGLEI KNGQALVALD KHWHLGCFKC KTCGKLLNAE YISKDGLPYC EADYHSKFGI RCDGCEKYIT GRVLEAGEKH YHPSCALCVR CGQMFSEGEE MYLQGSSIIWH PACRQAARTE DSKKETRTSS ESIVSVPASS TSGSPSRVIY AKLGDEILDY RDLAALPKNK AIYNIDRPDM ISYSPYISHS AVGDRQSYGE GDQDDRSYKQ CRTSSPSSAG SVSLGHYTPT SRSPQHYSRP GSESGRSTPS LSVHSDSRPP SSTYQQAPRH FHVPDTGVKD NIYRKPIIYK QHAARRLDVE DSSFQDSRK KTTWLLLKGD ADTRTNSPDL DSQSLSLSSG TDQEPLQRMA GDSLYSRFPY SKPDTLPGPR KDGLDLRAN LAPCGADPDA SWGTYREYKIY PYDSLIVTNR IRVKLPKDVD RTRLERHLSP EEFQEVEFGMS IEFDRLALW KRNDLKKKAL LF
Source	E.coli
Target Names	Ablim2
Protein Names	Recommended name: Actin-binding LIM protein 2 Short name= abLIM-2 Alternative name(s): Actin-binding LIM protein family member 2
Expression Region	1-612
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.