



Recombinant Mouse Eukaryotic translation initiation factor 1A, X-chromosomal (Eif1ax)

Product Code	CSB-BP807382MO
Abbreviation	Eif1ax
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BMJ3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPKKNKGKGGK NRRRGKNENE SEKRELVFKE DGQEYAQVIK MLGNRLEAM CFDGVKRLCH IRGKLRKKVW INTSDIILVG LRDYQDNKAD VILKYNAD EA RSLKAYGELP EHAKINETDT FGPGDDDEIQ FDDIGDDDED IDDI
Source	Baculovirus
Target Names	Eif1ax
Protein Names	Recommended name: Eukaryotic translation initiation factor 1A, X-chromosomal Short name= eIF-1A X isoform Alternative name(s): Eukaryotic translation initiation factor 4C Short name= eIF-4C
Expression Region	1-144
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein similar to eukaryotic translation initiation factor 1A (EIF1A). EIF1A is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNA ⁱ and eIF3) to the 5 end of capped RNA.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.