



Recombinant Mouse Protein NDNF (Ndnf)

Product Code	CSB-MP807450MO
Abbreviation	Ndnf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8C119
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	Q KLPTRDEELF QMQIRDKEFF HDSSVIPDGA EVSSYLFRDT PRRYFFMVVEE DNTPLSVTVT PCDAPLEWKL SLQELHEGSS ADGSGDPELL DQQKQQMTDV EGTELSYKKG NDVEYFLSSS SPSGLYQLEL LSTEKDTDFK VYATTTPESD QPYPELPYDP RVDVTSFGRT TVTLAWKPSP TASILKQPIE YCVVINKEHN FKSLCAAETK MNADDAFMVA PKPGLDFNPF DFAHFGFPTD NLGKDRSLLA KPSPKVGRHV YWRPKVDIQK ICIGNKNIFT VSDLKPDTQY YFDVFMVNTN TNMSTAYVGA FVRTKEEAKQ KTVELKDGRV TDVFKRKGG KFLRFAPVSS HQKVTFEHS CMDAVQVQVR RDGRLLLSQN VEGIRQFQLR GKPKGKYLIR LKGNRKGASK LKILATTRPS KHAFPSLPED TRIKAFDKLR TCSSVTVAWL GTQERRKFCI YRKEVDGNYS EDQKRREQNQ CLGPDTRKKS EKVLCYFHS QNLQKAVTTE TIRDLQPGKS YLLDVYVVGH GGHSVYQSK IVKTRKVC
Source	Mammalian cell
Target Names	Ndnf
Protein Names	Recommended name: Protein NDNF Alternative name(s): Epidermacan Neuron-derived neurotrophic factor
Expression Region	20-568
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.