



Recombinant Mouse Protein lin-9 homolog (Lin9)

Product Code	CSB-BP807470MO
Abbreviation	Lin9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8C735
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AELDQLPDE SSSAKALVSL KEGSLNTWN EKYSLSLQKTP VWKGRNAGPA VEMPFRNSKR SRLFSEDDDR QINTKSPKRN QRVAMIPQKF TATMSTPDKK ASQKIGFRLR NLLKLPKAHK WCIYEFYNSN IDKPLFEGDN DFCVCLKESF PNLKTRKLTR VEWGKIRRLM GKPRRCSSAF FEEERSALKQ KRQKIRLLQQ RKVADVSQFK DLPDEIPLPL VIGTKVTARL RGIHDGLFTG QIDAVDTLNA TYRVTFDRTG LGHTHTIPDYE VLSNEPHETM PISAFGQKQR PSRFFMTPPR LHYTPPLQSP ITDGDPLLQ SPWRSKVS GS DTETLGGFPV EFLIQVTKLS KILMIKKEHI KKLREMNT EA EKLKSYMPI GIEFQRRYAT IVLELEQLNK DLNKVLHKVQ QYCYELAPDQ GLQPADQPTD MRRRC EEEAQ EIVRQANSAS GQPCVENENL TDLISRLTAI LLQIKCLAEG GDLNSFEFKS LTDSLNDIKN TIDASNISCF QNNVEIHVAH IQSGLSQMGN LHAF AANNTN RD
Source	Baculovirus
Target Names	Lin9
Protein Names	Recommended name: Protein lin-9 homolog Short name= mLin-9 Alternative name(s): TUDOR gene similar 1 protein Type I interferon receptor beta chain-associated protein
Expression Region	2-542
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.