



Recombinant Mouse Cdc42-interacting protein 4 (Trip10)

Product Code	CSB-EP807544MO-B
Abbreviation	Trip10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CJ53
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MDWGTTELWDQ FEVLERHTQW GLDLLDKYVK FVKERAEVEQ AYAKQLRSLV KKYLPKRPTK DDPEVKFSQQ QSFVQLLQEV NDFAGQRELV AESLGIRVCL ELAKYSQEMK QERKMHFQEG RRAQQQLENG FKQLENSKRK FERDCREAEK AAHTAERLDQ DINATKADVE KAKQQAHLRN HMAEESKNEY AAQLQRFNRD QAHFYFSQMP QIFDKLQDMD ERRATRLGAG YGLLSEAELO VVPIIGKCLE GMKVAAESVD AKNDSQVLIE LHKSGFARPG DLEFEDFSQV INRVPSDSSL GTPDGRPELR AASSRSRAKR WPFGKKNKPR PPSLSLLGGH LPSTLSDGPS SPRSGRDPLA ILSEISKSVK PRLASFRSFR GGRGTVATED FSHLPPEQQR KRLQQQLEER NRELQKEEDQ REALKMKMDV YEKTPQMGGDP ASLEPRIAET LGNIERLKLE VQKYEAWLAE AESRVLSNRG DSLRSRHARPP DPPTTAPPDS SSSSTNSGSQ DNKESSEEP PSEGQDTPIY TEFDEDFEEP ASPIGQCVAI YHFEGSSEGT VSMSEGEDLS LMEEDKGDGW TRVRRKQGAE GYVPTSYLRLV TLN
Source	E.coli
Target Names	Trip10
Protein Names	Recommended name: Cdc42-interacting protein 4 Alternative name(s): Thyroid receptor-interacting protein 10 Short name= TR-interacting protein 10 Short name= TRIP-10
Expression Region	1-603
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.