



Recombinant Glutamate racemase (murl)

Product Code	CSB-YP808169EGX
Abbreviation	murl
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FB89
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MATKLQDGNI PCLAATPSEP RPTVLVFDSDG VGGLSVYDEI RHLLPDLHYI YAFDNVAFPY GEKSEVFIVE RVVEIVTAVQ ERYPLALAVV ACNTASTVSL PALREKDFDP VVGVVPAIKP AARLTANGIV GLLATRGTVK RSYTHELIAR FANECQIEML GSAEMVELAE AKLHGEDVSL DALKRILRPW LRMKEPPDTV VLGCTHFPLL QEELLQVLPE GTRLVDSGAA IARRTAWLLE HEAPDAKSAD ANIAFCMAMT PEAEQLLPVL QRYGFETLEK LAVLD
Source	Yeast
Target Names	murl
Protein Names	Recommended name: Glutamate racemase EC= 5.1.1.3
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>