



# Recombinant Ketol-acid reductoisomerase (ilvC)

<b>Product Code</b>	CSB-BP808174EGX
<b>Abbreviation</b>	ilvC
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q8FBR2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ANYFNTLNL RQQLAQLGKC RFMGRDEFAD GASYLQGKKV VIVGCGAQGL NQGLNMRDSG LDISYALRKE AIAEKRASWR KATENGFKVG TYEELIPQAD LVVNLTPDKQ HSDVVRTVQP LMKDGAALGY SHGFNIVEVG EQIRKDITVV MVAPKCPGTE VREEYKRGFG VPTLIAVHPE NDPKGEGMAI AKAWAAATGG HRAGVLESSF VAEVKSDLMG EQTILCGMLQ AGSLLCFDKL VEEGTPAYA EKLIQFGWET ITEALKQGGI TLMMDRLSNP AKLRAYALSE QLKEIMAPLF QKHMDDIISG EFSSGMMADW ANDDKLLTW REETGKTAFE TAPQYEGKIG EQEYFDKGV LMIAMVKAGVE LAFETMVDSG IIEESAYYES LHELPLIANT IARKRLYEMN VVISDTAEYG NYLFSYACVP LLKPFMAELQ PGDLGKAIFE GAVDNAQLRD VNEAIRSHAI EQVGKKLRGY MTDMKRIAVA G
<b>Source</b>	Baculovirus
<b>Target Names</b>	ilvC
<b>Protein Names</b>	Recommended name: Ketol-acid reductoisomerase EC= 1.1.1.86 Alternative name(s): Acetohydroxy-acid isomeroeductase Alpha-keto-beta-hydroxylacil reductoisomerase
<b>Expression Region</b>	2-491
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>