



Recombinant Xylose import ATP-binding protein XylG (xylG)

Product Code	CSB-MP808178EGX
Abbreviation	xylG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FCE2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MPYLLEMKNI TKTFGSVKAI DNVSLRLNAG EIVSLCGENG SGKSTLMKVL CGIYPHGSYE GEIIFAGEEI QASHIRDTER KGIAIIHQEL ALVKELTVLE NIFLGNETH NGIMDYDLMT LRCQKLLAQV SLSISPDTRV GDLGLGQQQL VEIAKALNKQ VRLILDEPT ASLTEQETSV LLDIIRDLQQ HGIACIYISH KLNEVKAISD TICVIRDGQH IGTRDAAGMS EDDIITMMVG RELTALYPNE PHTTGDEILR IEHLTAWHPV NRHIKRVNDV SFSLKRGEIL GIAGLVGAGR TETIQCLFGV WPGQWEGKIY IDGKQVDIRN CQQAIAQGIA MVPEDRKRDRG IVPMAVGKN ITLAALKKFT GSISQLDDAA EQKCILESIQ QLKVKTSPPD LAIGRLSGGN QKAILARCL LLNPRILILD EPTRGIDIGA KYEIKLINQ LVQQGIAVIV ISSELPEVLG LSDRVLVMHE GKLKANLINH NLTQEQVMEA ALRSEHHVEK QSV
Source	Mammalian cell
Target Names	xylG
Protein Names	Recommended name: Xylose import ATP-binding protein XylG EC= 3.6.3.17
Expression Region	1-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.