



Recombinant Uncharacterized protein ygeW (ygeW)

Product Code	CSB-EP808189EGX
Abbreviation	ygeW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FE91
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MMKTVNELIK DINSLTSHLH EKDFLLTWEQ TPDELKQVLD VAAALKALRA ENISTKVFNS GLGISVFRDN STRTRFSYAS ALNLLGLAQQ DLDEGKSQIA HGETVRETAN MISFCADAIG IRDDMYLGAG NAYMREVGAA LDDGYKQGV PQRPALVNLQ CDIDHPTQSM ADLAWLREHF GSLENLKGKK IAMTWAYSPS YGKPLSVPQG IIGLMTRFGM DVTLAHPEGY DLIPDVVEVA KNNAKASGGS FRQVTSMEEA FKDADIVYPK SWAPYKVMEE RTELLRANDH EGLKALEKQC LAQNAQHKDW HCTEEMMELT RDGEALYMHCLPADISGVSC KEGEVTEGVF EKYRIATYKE ASWKPYIIAA MILSRKYAKP GALLEQLLKE AQERVK
Source	E.coli
Target Names	ygeW
Protein Names	Recommended name: Uncharacterized protein ygeW
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.