



Recombinant Phosphoadenosine phosphosulfate reductase (cysH)

Product Code	CSB-EP808191EGX-B
Abbreviation	cysH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FEI9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	SKLDLNALN ELPKVDRILA LAETNAQLEK LDAEGRVAWA LDNLPGEYVL SSSFGIQA AV SLHLVNQIRP DIPVILTD TG YLFPETYRFI DELTDK LKLN LKVYRATESA AWQEARYGKL WEQGVEGIEK YNDINKVEPM NRALKELNAQ TWFAGLRREQ SGS RANLPVL AIQRGVFKVL PIIDWDNRTI YQYLQKHGLK YHPLWDEGYL SVGDTH TTRK WEPGMAEEET RFFGLKRECG LHEG
Source	E.coli
Target Names	cysH
Protein Names	Recommended name: Phosphoadenosine phosphosulfate reductase EC=1.8.4.8 Alternative name(s): 3'-phosphoadenylylsulfate reductase PAPS reductase, thioredoxin dependent PAPS sulfotransferase PAdoPS reductase
Expression Region	2-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.