



Recombinant S-ribosylhomocysteine lyase (luxS)

Product Code	CSB-EP808193EGX
Abbreviation	luxS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FEP8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	PL LDSFTVD HTRMEAPAVR VAKTMNTPHG DAITVFDLRF CVPNKEVMPE RGIHTLEHLF AGFMRNHLNG NGVEIIDISP MGCRTGFYMS LIGTPDEQRV ADAWKAAMED VLKVQDQIQI PELNVYQCGT YQMHSLQEAQ DIARNILERD VRINSNEELA LPKEKLQELH I
Source	E.coli
Target Names	luxS
Protein Names	Recommended name: S-ribosylhomocysteine lyase EC= 4.4.1.21 Alternative name(s): AI-2 synthesis protein Autoinducer-2 production protein LuxS
Expression Region	2-171
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.