



Recombinant tRNA sulfurtransferase (thil)

Product Code	CSB-BP808236EGX
Abbreviation	thil
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FKB7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MKFIKLFPE ITIKSQSVRL RFIKILTGNI RNVLKHYDET LAVVRHWDNI EVRAKDENQR LAIRDALTRI PGIHHILEVE DVPFTDMHDI FEKALVQYRD QLEGKTFCVR VKRRGKHDFS SIDVERYVGG GLNQHIESAR VKLTNPEVTV HLEVEDDRLL LIKGRYEGIG GFPIGTQEDV LSLISGGFDS GVSSYMLMRR GCRVHYCFFN LGGAAHEIGV RQVAHYLWNR FGSSHRVRFV AINFEPVVGGE ILEKIDDGQM GVILKRMMVR AASKVAERYG VQALVTGEAL GQVSSQTLTN LRLIDNVSDT LILRPLISYD KEHIINLARQ IGTEDFARTM PEYCGVISKS PTVKAVKSKI EAEEEKDFDS ILDKVVEEAN NVDIREIAQQ TEQEVVEVET VNGFGPNDVI LDIRSIDEQE DKPLKVEGID VVSLPFYKLS TKFGDLDQSK TWLLWCERGV MSRLQALYLR EQGFNNVKVY RP
Source	Baculovirus
Target Names	thil
Protein Names	Recommended name: tRNA sulfurtransferase EC= 2.8.1.4 Alternative name(s): Sulfur carrier protein ThiS sulfurtransferase Thiamine biosynthesis protein Thil tRNA 4-thiouridine synthase
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.