



# Recombinant *Oryza sativa* subsp. japonica Putative zinc finger CCCH domain-containing protein 48 (Os07g0139000, LOC\_Os07g04650)

<b>Product Code</b>	CSB-EP808410OFG-B
<b>Abbreviation</b>	Os07g0139000, LOC_Os07g04650
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8GVZ8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADSEDGDHV ATTTTEQRYD DDGHLVPSSG GQEEEGSGGR DVVVPGGHVA EDYRSGVGVP VGRDAGGATS SPPQPVHVTP SILVGSIHAP VFQGELVGMK FGVGSGSMGA GTSATRRLPA TGFGALPTSS MAEDSADHAD DDHLAEEEEEE EEEHYIDDG PLVPSSGGQE EEGSGGRHVF VPGGHDGEED HPDDLVDLD LDLLVDGVVG PVPGGHLNAD APAFVPTTRG RQDLYSALSS SAPAAGYRYR HYITSSALAE AGHVSPFLGL PYATAFDSPL DRELVGPSSA PPPCSAASRA WLVRCSSPLS DSEWTRRSIL AREAHTPAS TVTGRGRFEF VPIPGAPYAP PPSFAPIAAG AGPAARPLQQ LAFGLEEHKT KLCAEYYSRG LGCPRGNTCK YAHGEDDLRL VVAVSSLADA GEGSSSSDSS FAALGGEDKY KTKLCKTFTS GGLCLFAANC RFAHGEVELG KKEPCWYFFS GQTCPRGDTC GFRHSY
<b>Source</b>	E.coli
<b>Target Names</b>	Os07g0139000
<b>Protein Names</b>	Recommended name: Putative zinc finger CCCH domain-containing protein 48 Short name= OsC3H48
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.