



Recombinant *Oryza sativa* subsp. japonica Putative zinc finger CCCH domain-containing protein 48 (Os07g0139000, LOC_Os07g04650)

Product Code	CSB-MP8084100FG
Abbreviation	Os07g0139000, LOC_Os07g04650
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GVZ8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MADSEDGDHV ATTTTEQRYD DDGHLVPSSG GQEEEGSGGR DVVVPGGHVA EDYRSGVGVP VGRDAGGATS SPPQPVHVTP SILVGSIHAP VFQGELVGMK FGVGSGSMGA GTSATRRLPA TGFGALPTSS MAEDSADHAD DDHLAEEEEEE EEEHYIDDG PLVPSSGGQE EEGSGGRHVF VPGGHDGEED HPDDLVDLD LDLLVDGVVG PVPGGHLNAD APAFVPTTRG RQDLYSALSS SAPAAGYRYR HYITSSALAE AGHVSPFLGL PYATAFDSPL DRELVGPSSA PPPCSAASRA WLVRCSSPLS DSEWTRRSIL AREAHTPAS TVTGRGRFEF VPIPGAPYAP PPSFAPIAAG AGPAARPLQQ LAFGLEEHKT KLCAEYYSRG LGCPRGNTCK YAHGEDDLRL VVAVSSLADA GEGSSSSDSS FAALGGEDKY KTKLCKTFTS GGLCLFAANC RFAHGEVELG KKEPCWYFFS GQTCPRGDTC GFRHSY
Source	Mammalian cell
Target Names	Os07g0139000
Protein Names	Recommended name: Putative zinc finger CCCH domain-containing protein 48 Short name= OsC3H48
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.