



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein SKIP25 (SKIP25)

Product Code	CSB-YP808415DOA
Abbreviation	SKIP25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GX29
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEKCLKRRES MSTTAAESPP AKRRRTVTGN ENSALIEGLP DHISEICLSL VHRPSLLSAV CTRWRRLLYS PEFPSFPSLY ALFVDSTSDT GRVNPSVRFM CFNPVSSKWKY PLPPPPDPP LHRILYRHPS FISFNLPIQC VSAAGKLIL AGSNQQLSPA ISHPLIFDPI SSSWSSGPRI GSPRRWCATG ACDGAIYIAS GISSQFSSTV AKSVEKLDLT EQNRNNHRFN WEKLRDMRD L RFSREIDA V GYRRKLLMVN VKGDAVKEGA IYDVVKDDWE PMPEEMLVGW RGPVAAMEEE ILYSVDERRG TVRKYDDEKR EWREVVVVEG GEEMKLGATQ VTADSGKLCV VTGDGKIVVV DVAAEPAKIW NVEIPDGLEP VSVHVLPRMS QPNFC
Source	Yeast
Target Names	SKIP25
Protein Names	Recommended name: F-box/kelch-repeat protein SKIP25 Alternative name(s): SKP1-interacting partner 25
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.