



Recombinant Arabidopsis thaliana Transcription factor bHLH63 (BHLH63)

Product Code	CSB-EP808425DOA-B
Abbreviation	BHLH63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GY61
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MNGAIGGDLL LNFPMDSVLE RQRAHLKYLN PTFDSPLAGF FADSSMITGG EMDSYLSTAG LNLPMMYGET TVEGDSRLSI SPETTLGTGN FKRRKFDтет KDCNEKKKKM TMNRDDLVEE GEEEEKSKITE QNNGSTKSIK KMKHKAKKEE NNFSNDSSKV TKELEKTDYI HVRARRGQAT DSHSIAERVR REKISERMKF LQDLVPGCDK ITGKAGMLDE IINYVQSLQR QIEFLSMKLA IVNPRPDFDM DDIFAKEVAS TPMTVVPSPE MVLSGYSHEM VHSGYSSEMV NSGYLHVNPM QQVNTSSDPL SCFNNGEAPS MWDSHVQONLY GNLGV
Source	E.coli
Target Names	BHLH63
Protein Names	Recommended name: Transcription factor bHLH63 Alternative name(s): Basic helix-loop-helix protein 63 Short name= AtbHLH63 Short name= bHLH 63 Transcription factor EN 84 bHLH transcription factor bHLH063
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.