



Recombinant Arabidopsis thaliana KH domain-containing protein At3g08620 (At3g08620)

Product Code	CSB-MP808430DOA
Abbreviation	At3g08620
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GYR4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSGLYNYNNF SPSRAASPQI RTPSSDVDSQ YISQLLAEHQ KLGPFMQVLP ICSRLLNQEI FRITGMMPNQ GFTDFDRLRH RSPSPMASPN LMSNVSGGGL GGWNGLPPE IGPHGMAME WQGAPASPSS YPVKRILRLD LPVDTPNFN FVGRLLGPRG NSLKRVEATT GCRVYIRGKG SIKDPEKEEK LKGKPGYEHL NEQLHILIEA DLPIDIVDIK LRQAQEIIEE LVKPVDESQD YIKRQQLREL ALLNSNLREN SPGPSGSVSP FNSNAMKRPK TGR
Source	Mammalian cell
Target Names	At3g08620
Protein Names	Recommended name: KH domain-containing protein At3g08620 Alternative name(s): Quaking-like protein 4
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.