



Recombinant Arabidopsis thaliana Transcription factor bHLH128 (BHLH128)

Product Code	CSB-MP808437DOA
Abbreviation	BHLH128
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H102
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MYQSSSSTSS SSQRSSLPGG GGLIRYGSAP GSFLNSVVDE VIGGGSSNAR DFTGYQPSSD NFIGNFFTGA ADSSSLRSDS TTCGVNNSSD GQKQLGNNNN NNSNKDIFLD RSYGGFNEIS QQHKSNDIGG GNSSGSYSLA RQRSSPADFF TYLASDKNNF SLNQPTSDYS PQGGSNGGRG HSRLKSQLSF TNHDSLARIN EVNETPVHDG SGHSFSAASF GAATTDSDWDD GSGSIGFTVT RPSKRSKDMD SGLFSQYSLP SDTSMNYMDN FMQLPEDSVP CKIRAKRGCA THPRSAIERE RRTRISGKLG KLQDLVPNMD KQTSYSDMLD LAVQHIKGLQ HQLQNLKKDQ ENCTCGCSEK PS
Source	Mammalian cell
Target Names	BHLH128
Protein Names	Recommended name: Transcription factor bHLH128 Alternative name(s): Basic helix-loop-helix protein 128 Short name= AtbHLH128 Short name= bHLH 128 Transcription factor EN 74 bHLH transcription factor bHLH128
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.