



Recombinant Dog Glutamine synthetase (GLUL)

Product Code	CSB-BP808499DO
Abbreviation	GLUL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8HZM5
Product Type	Recombinant Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	>85% (SDS-PAGE)
Sequence	ATSASSHLN KGIKQVYMSL PQGEKVQAMY IWIDGTGEGL RCKTRTL DSE PKGVEELPEW NFDGSSTFQS EGSNSDMYLV PAAMFRDPFR KDPNKLVFCE VFKYNRKP AE TNL RHTCKRI MDMVSNQHPW FGMEQEY TLM GTDGHPFGWP SNGFPGPQGP YYCGVGADKA YGRDIVEAHY RACLYAGIKI AGTNAEVMPA QWEFQIGPCE GIDMGDHLWV ARFILHRVCE DFGVIATFDP KPIPGNWNGA GCHTNFSTKA MREENGLKYI EESIEKLSKR HQYHIRAYDP KGGLDNARRL TGFHETSNIN DFSAGVANRG ASIRIPRTVG QEKKGYFEDR RPSANCDPFS VTEALIRTCL LNETGDEPFQ YKN
Source	Baculovirus
Target Names	GLUL
Protein Names	Recommended name: Glutamine synthetase Short name= GS EC= 6.3.1.2 Alternative name(s): Glutamate decarboxylase EC= 4.1.1.15 Glutamate-- ammonia ligase
Expression Region	2-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.