



# Recombinant Human Protein hid-1 homolog (C17orf28)

<b>Product Code</b>	CSB-BP808538HU
<b>Abbreviation</b>	C17orf28
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IV36
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GSTD SKLNF RKAVIQLTTK TQPVEATDDA FWDQFWADTA TSVQDVFALV PAAEIRAVRE ESPSNLATLC YKAVEKLVQG AESGCHSEKE KQIVLNC SRL LTRVLPYIFE DPDWRGFFWS TVPGAGRGGQ GEEDDEHARP LAESLLAIA DLLFCPDFTV QSHRRSTVDS AEDVHSLDSC EYIWEAGVGF AHSPQPNYIH DMNRMELLKL LLTCFSEAMY LPPAPESGST NPWVQFFCST ENRHALPLFT SLLNTVCAYD PVGYGIPYNH LLFSDYREPL VEEAAQVLIV TLDHDSASSA SPTVDGTTTG TAMDDADPPG PENLNFVNYLS RIHREEDFQF ILKGIARLLS NPLLQTYLPN STKKIQFHQE LLVLFWKLCD FNKKFLFFVL KSSDVL DILV PILFFLNDAR ADQSRVGLMH IGVFILLLS GERNFGVRLN KPYSIRVPMD IPVFTGTHAD LLIVVFHKII TSGHQRLQPL FDCLLTIVVN VSPYLKSLSM VTANKLLHLL EAFSTTWFLF SAAQNHHLVF FLLEVFNNII QYQFDGNSNL VYAIIRKRSI FHQLANLPTD PPTIHKALQR RRRTPPELSR TGSQEGTSME GSRPAAPAEP GTLKTSLVAT PGIDKLTEKS QVSEDGTLRS LEPEPQQSLE DGSPAKGEPS QAWREQRRPS TSSASGQWSP TPEWVLSWKS KLPLQTIMRL LQVLVPQVEK ICIDKGLTDE SEILRFLQHG TLVGLLPVPH PILIRKYQAN SGTAMWFRTY MWGVIYLRNV DPPVWYDTDV KLFEIQRV
<b>Source</b>	Baculovirus
<b>Target Names</b>	HID1
<b>Protein Names</b>	Recommended name: Protein hid-1 homolog Alternative name(s): Down-regulated in multiple cancers 1
<b>Expression Region</b>	2-788
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.