



Recombinant Human Protein hid-1 homolog (C17orf28)

Product Code	CSB-YP808538HU
Abbreviation	C17orf28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IV36
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GSTD SKLNF RKAVIQLTTK TQPVEATDDA FWDQFWADTA TSVQDVFALV PAAEIRAVRE ESPSNLATLC YKAVEKLVQG AESGCHSEKE KQIVLNC SRL LTRVLPYIFE DPDWRGFFWS TVPGAGRGGQ GEEDDEHARP LAESLLAIA DLLFCPDFTV QSHRRSTVDS AEDVHSLDSC EYIWEAGVGF AHSPQPNYIH DMNRMELLKL LLTCFSEAMY LPPAPESGST NPWVQFFCST ENRHALPLFT SLLNTVCAYD PVGYGIPYNH LLFSDYREPL VEEAAQVLIV TLDHDSASSA SPTVDGTTTG TAMDDADPPG PENLNFVNYLS RIHREEDFQF ILKGIARLLS NPLLQTYLPN STKKIQFHQE LLVLFWKLCD FNKKFLFFVL KSSDVL DILV PILFFLNDAR ADQSRVGLMH IGVFILLLS GERNFGVRLN KPYSIRVPMD IPVFTGTHAD LLIVVFHKII TSGHQRLQPL FDCLLTIVVN VSPYLKSLSM VTANKLLHLL EAFSTTWFLF SAAQNHHLVF FLLEVFNNII QYQFDGNSNL VYAIIRKRSI FHQLANLPTD PPTIHKALQR RRRTPPELSR TGSQEGTSME GSRPAAPAEP GTLKTSLVAT PGIDKLTEKS QVSEDGTLRS LEPEPQQSLE DGSPAKGEPS QAWREQRRPS TSSASGQWSP TPEWVLSWKS KLPLQTIMRL LQVLVPQVEK ICIDKGLTDE SEILRFLQHG TLVGLLPVPH PILIRKYQAN SGTAMWFRTY MWGVIYLRNV DPPVWYDTDV KLFEIQRV
Source	Yeast
Target Names	HID1
Protein Names	Recommended name: Protein hid-1 homolog Alternative name(s): Down-regulated in multiple cancers 1
Expression Region	2-788
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.