



Recombinant Human Tryptophan 5-hydroxylase 2 (TPH2)

Product Code	CSB-MP809000HU
Abbreviation	TPH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IWU9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQPAMMMFSS KYWARRGFSL DSAVPEEHQL LGSSTLNKPN SGKNDDKGNK GSSKREAATE SGKTAVVFSL KNEVGGLVKA LRLFQEKRVN MVHIESRKSR RRSSEVEIFV DCECGKTEFN ELIQLLKFT TIVTLNPPEN IWTEEELEED VPWFPRKISE LDKCSHRVLM YGSELDADHP GFKDNVYRQR RKYFVDVAMG YKYGQPIPRV EYTEEETKTW GVVFRELSKL YPTHACREYL KNFPLLTKYC GYREDNVPQL EDVSMFLKER SGFTVRPVAG YLSPRDFLAG LAYRVFHCTQ YIRHGSDPLY TPEPDTCHEL LGHVPLLADP KFAQFSQEIG LASLGASDED VQKLATCYFF TIEFGLCKQE GQLRAYGAGL LSSIGELKHA LSDKACVKAF DPKTTCLQEC LITTFQEAYF VSESFEEAKE KMRDFAKSIT RPFVSVYFNPY TQSIEILKDT RSIENVVQDL RSDLNTVCDA LNKMNQYLGI
Source	Mammalian cell
Target Names	TPH2
Protein Names	Recommended name: Tryptophan 5-hydroxylase 2 EC= 1.14.16.4 Alternative name(s): Neuronal tryptophan hydroxylase Tryptophan 5-monoxygenase 2
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the pterin-dependent aromatic acid hydroxylase family. The encoded protein catalyzes the first and rate limiting step in the biosynthesis of serotonin, an important hormone and neurotransmitter. The human genome contains two related tryptophan hydroxylases, one on chromosome 11p15-p14 and one on chromosome 12q21. This gene is expressed predominantly in the brain stem. Mutations in this gene may be associated with psychiatric diseases such as bipolar affective disorder and major depression.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.