



Recombinant Human Uncharacterized protein C11orf53 (C11orf53)

Product Code	CSB-YP809008HU
Abbreviation	C11orf53
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXP5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPGSPVTSGY YGVRRSFLSD SDFHNSKQFS NDVYTSSVGK PFPCESSAGQ SHAALLEPYF PQEPYGDYRP PALTPNAGSL FSASPLPPLL PPPFPGDPAH FLFRDSWEQT LPDGLSQPDP VSADALLTLP PSTSCLSQLE SGSIAQHRGS SWGSSLAGAQ SYSLHALEDL HHTPGYPTPP PYPFTPMTV SNDLPPKVGP LSPDEEADTG SLHDPSPWVK EDGSIWGSY ECRRAY
Source	Yeast
Target Names	C11orf53
Protein Names	Recommended name: Uncharacterized protein C11orf53
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.