



Recombinant Human Protein AMN1 homolog (AMN1)

Product Code	CSB-EP809013HU
Abbreviation	AMN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IY45
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPRRRVSQL LDLCLWCFMK NISRYLTDIK PLPPNIKDRL IKIMSMQQGI TDSNISEILH PEVQTLDLRS CDISDAALLH LSNCRKLKKL NLNASKGNRV SVTSEGIKAV ASSCSYLHEA SLKRCCNLTD EGVVALALNC QLLKIIDLGG CLSITDVSLH ALGKNCPFLQ CVDFSATQVS DSGVIALVSG PCAKKLEEIH MGHCVNLTDG AVEAVLTYCP QIRILLFHGC PLITDHSREV LEQLVGPNKL KQVTWTVY
Source	E.coli
Target Names	AMN1
Protein Names	Recommended name: Protein AMN1 homolog
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.