



Recombinant Human Cytoskeleton-associated protein 2-like (CKAP2L)

Product Code	CSB-MP809015HU
Abbreviation	CKAP2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYA6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVGGPGPTAAA AVEERQRKLQ EYLAAGKGLK SQNTKPYLKS KNNCQNQPPS KSTIRPKNDV TNHVVLPVKP KRSISIKLQP RPPNTAGSQQ PKLEPPKLLG KRLTSECVSS NPYSKPSSKS FQQCEAGSST TGELSRKPVG SLNIEQLKTT KQQLTDQGNG KCIDFMNNIH VENESLDNFL KETNKENLLD ILTEPERKPD PKLYTRSKPK TDSYNQTKNS LVPKQALGKS SVNSAVLKDR VNKQFVGETQ SRTFPVKSQQ LSRGADLARP GVKPSRTVPS HFIRTLSKVQ SSKKPVVKNI KDIKVNRSQY ERPNETKIRS YPVTEQRVKH TKPRTYPSLL QGEYNNRHPN IKQDQKSSQV CIPQTSCVLQ KSKAISQRPN LTVGRFNSAI PSTPSIRPNG TSGNKHNNNG FQQKAQTLDS KLKKAQPQNH FLNKTAPKTQ ADVTTVNGTQ TNPNIKKKAT AEDRRKQLEE WQKSKGKTYK RPPMELKTKR KVIKEMNISF WKSIEKEEEE KKAQLELSSK INNTLTECLN LIEGGVPSNE ILNILSSIPE AEKFAKFWIC KAKLLASKGT FDVIGLYEEA IKNGATPIQE LRKVVLNILQ DSNRTTEGIT SDSLVAETSI TSVEELAKKM ESVKSCLSPK EREQVTATPR IAKAEQHNYP GIKLQIGPIP RINGMPEVQD MKFITPVRRS SRIERAVSRY PEMLQEHDLV VASLDELLEV EETKCFIFRR NEALPVTLGF QTPES
Source	Mammalian cell
Target Names	CKAP2L
Protein Names	Recommended name: Cytoskeleton-associated protein 2-like
Expression Region	1-745
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.