



# Recombinant Human UPF0704 protein C6orf165 (C6orf165)

<b>Product Code</b>	CSB-BP809020HU
<b>Abbreviation</b>	C6orf165
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYR0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPPTQAESVI RSIIREIGQE CAAHGEIVSE TLIAFMVKAV VLDPSNGFNM DRTLKMSDVQ NLVKLCMTRL LDTKNPSLDT IKMQVYFDMN YTNRVEFLEE HHRVLESRLG SVTREITDNR ACAKEELES YRKIISYVLL RSGLSPTDI KTVREVTAAL QSVFPQAE LG TFLTLSKKDK ERQLKELTMI VTGIRLFNRD CGKGGEGIDD LPAVLHVAIP ATMQHIDYQL ETARSQVYRY TAILEKAAND PLMRAELQPY MLKEALYNIR QYEVFLQIIL SDIITGAQEV EMMTKQLGAH LEQLKMTIKS KIAVPTSQVF PIFIALSTLW TSLQDETIVV GVLSNLFTHI QPFLGAHELY FPERVMQCHL NGATVKTDVC RMKEHMEDRV NVADFRKLEW LFPETTANFD KLLIQYRGFC AYTFAATDGL LLPGNPAIGI LKYKEKYYTF NSKDAAYSFA ENPEHYIDIV REKAKKNTTEL IQLLELHQQF ETFIPYSQMR DADKHYIKPI TKCESSTQTN THILPPTIVR SYEWNEWELR RKAIKLANLR QKVTHSVQTD LSHLRRENC S QVYPPKDTST QSMREDSTGV PRPQIYLAGL RGGKSEITDE VKVNLTRDVD ET
<b>Source</b>	Baculovirus
<b>Target Names</b>	CFAP206
<b>Protein Names</b>	Recommended name: UPF0704 protein C6orf165
<b>Expression Region</b>	1-622
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.