



# Recombinant Human Katanin p60 ATPase-containing subunit A-like 2 (KATNAL2)

<b>Product Code</b>	CSB-YP809022HU
<b>Abbreviation</b>	KATNAL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYT4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MELSYQTLKF THQAREACEM RTEARRKNLL ILISHYLTQE GYIDTANALE QETKLGLRRF EVCDNIDLET ILMEYESYYF VKFQKYPKIV KKSSDTAENN LPQRSRGKTR RMMNDSCQNL PKINQQRPRS KTTAGKTGDT KSLNKEHPNQ EVVDNTRLES ANFGLHISRI RKDSGEENAH PRRGQIIDFQ GLLTDAIKGA TSELALNTFD HNPDPSEKLL KPLSAFIGMN SEMRELAADV SRDIYLNPN IKWNDIIGLD AAKQLVKEAV VYPIRYPQLF TGILSPWKGL LLYGPPGTGK TLLAKAVATE CKTTFNISA STIVSKWRGD SEKLVKRVLFE LARYHAPSTI FLDELESVMS QRGTASGGEH EGSLRMKTEL LVQMDGLARS EDLVFVLAAS NLPWELDCAM LRRLEKRILV DLPSREARQA MIYHWLPPVS KSRALELHTE LEYSVLSQET EGYSGSDIKL VCREAAMPV RKIFDALENH QSESSDLPRI QLDIVTTADF LDVLTHTKPS AKNLAQRYSD WQREFESV
<b>Source</b>	Yeast
<b>Target Names</b>	KATNAL2
<b>Protein Names</b>	Recommended name: Katanin p60 ATPase-containing subunit A-like 2 Short name= Katanin p60 subunit A-like 2 EC= 3.6.4.3 Alternative name(s): p60 katanin-like 2
<b>Expression Region</b>	1-538
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.