



Recombinant Human Dead end protein homolog 1 (DND1)

Product Code	CSB-BP809023HU
Abbreviation	DND1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYX4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQSKRDCELW CERVNPENKA ALEAWVRETG IRLVQVNGQR KYGGPPPGWV GSPPPAGSEV FIGRLPQDVY EHQLIPLFQR VGRLYEFRLM MTFSGLNRRGF AYARYSSRRG AQAAIATLHN HPLRPSCPLL VCRSTEKCEL SVDGLPPNLT RSALLLALQP LGPGLQEARL LPSPGPAPGQ IALLKFSSHR AAAMAKKALV EGQSHLCGEQ VAVEWLKPD LKQRLRQQLVG PFLRSPQPEG SQLALARDKL GFQGARATLQ LLCQRMKLG S PVFLTCLGI GPAGWHRFWY QVVIPGHPVP FSGLIWVLT LDGRDGHEVA KDAVSVRL LQ ALSESGANLL WSAGAEAGTM VKQ
Source	Baculovirus
Target Names	DND1
Protein Names	Recommended name: Dead end protein homolog 1 Alternative name(s): RNA-binding motif, single-stranded-interacting protein 4
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.