



Recombinant Human Protein-associating with the carboxyl-terminal domain of ezrin (SCYL3)

Product Code	CSB-MP809026HU
Abbreviation	SCYL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZE3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GSENSALKS YTLREPPFTL PSGLAVYPAV LQDGKFASVF VYKRENEDKV NKAAKHLKTL RHPCLLRFLS CTVEADGIHL VTERVQPLEV ALETLSAEV CAGIYDILLA LIFLHDRGHL THNNVCLSSV FVSEDGHWKL GGMETVCKVS QATPEFLRSI QSIRDPASIP PEEMSPEFTT LPECHGHARD AFSFGTLVES LLTILNEQVS ADVLSSFQQT LHSTLLNPIP KCRPALCTLL SHDFFRNDFL EVVNFLKSLT LKSEEEKTEF FKFLLDREVSC LSEELIASRL VPLLLNQLVF AEPVAVKSFL PYLLGPKKDH AQGETPCLLS PALFQSRVIP VLLQLFEVHE EHVRMVLLSH IEAYVEHFTQ EQLKKVILPQ VLLGLRDTSD SIVAITLHSL AVLVSLLGPE VVVGGERTKI FKRTAPSFTK NTDLSLEDSP MCVVCSHHSQ ISPILENPFS SIFPKCFFSG STPINSKKHI QRDYYNTLLQ TGDPFSSQPIK FPINGLSDVK NTSEDSSENF SSSKKSEEW P DWSEPEEPEN QTVNIQIWPR EPCDDVKSQC TTLDVEESSW DDCEPSSLDT KVNPGGGITA TKPVTSGEQK PIPALLSLTE ESMPWKSSLP QKISLVQRGD DADQIEPPKV SSQERPLKVP SELGLGEEFT IQVKKKPVKD PEMDWFADMI PEIKPSAAFL ILPELRTEMV PKKDDVSPVM QFSSKFAAAE ITEGAEAGWE EEGELNWEDN NW
Source	Mammalian cell
Target Names	SCYL3
Protein Names	Recommended name: Protein-associating with the carboxyl-terminal domain of ezrin Alternative name(s): Ezrin-binding protein PACE-1 SCY1-like protein 3
Expression Region	2-742
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.